

The 1x1 tDCS device packages the most advanced stimulation technology in an intuitive user interface.

SIMPLE AND SMART

Set stimulation settings in seconds and prevent errors with simple to set and read controls.

DURATION

Current duration of up to 40 minutes of continuous stimulation.

INTENSITY

Easy adjustment of stimulation intensity to up to 5 mA. Applies a fixed ramp slope current output.




Whether you want to apply standard tDCS protocols or establish new ones, our scientific team is ready to work with you to determine optimal dose parameters and accessories.



contact@soterixmedical.com



www.soterixmedical.com



1x1 tDCS

INDUSTRY STANDARD

DEVICE FEATURES

TIME REMAINING

Displays time remaining in session.

SMARTscan

Continuous display of contact conditions prior, during, and after stimulation.

PRE-STIM TICKLE

Condition skin resistance prior to stimulation as needed.

TRUE CURRENT

Continuous display of actual current value being delivered regardless of device settings.

RELAX

Titrate output current in increments of 10uA.

ABORT

Terminate stimulation at any time. Current will ramp down to zero.

